

SteriSpine™ PS

Pedicle Screw System

Full traceability

Instruments and implants are fully traceable, creating a safer process for patients and reducing risk for hospitals and surgeons.

Eliminate Pre-Op Costs

Pre-op handling, wrapping and sterilization is eliminated, creating a more cost-effective process for the hospital.

Surgical Technique Flexibility

Universal instrument kit provides everything necessary for open and MIS surgical techniques

Always Surgery Ready

Instruments and implants are always clean, sterile, in perfect calibration and surgery-ready at a moment's notice, increasing OR staff efficiency

Procedural Efficiency

Preloaded screw and set screw assemblies create more simple and efficient OR experience.



Eliminate Cross-Contamination Risk

Single-use and sterile instrument kit eliminates the risk of instrument cross-contamination.

Eliminate Post-Op Costs

Post-op handling, cleaning and decontamination is eliminated, creating a more cost-effective process for the hospital

Always New

Sealed, sterile kits are always complete with never used before instrumentation, eliminating the risk of damaged, unclean or missing instruments.

Fast OR Turnaround

Lightweight, sterile and traceable kits enhance OR team efficiency by eliminating replenishment and processing delays between cases

Contact Safe Orthopaedics

For more information about our products and company, please visit www.safeortho.com or contact us today.

Western Europe / International

Safe Orthopaedics, SAS
Parc des Bellevues
Allée R. Luxembourg, Le Californie
95610 Eragny sur Oise
FRANCE
Ph: +33 1 34 21 12 43
Fax: +33 1 34 21 12 00

North America

Safe Orthopaedics, LLC
3175 Lenox Park Blvd.
Suite 108
Memphis, TN 38115
USA
Toll Free Ph: (855) 837-3959
Toll Free Fax: (855) 837-3958

SteriSpine™ PS

Pedicle Screw System



Instrumentation kit for Open or MIS Surgery

Ref : KITANCMIS

PSANC018
Guide Wire
x6



PSANC009
Awl
x1



PSANC010
Probe
x1




PSANC012
Dual Cannulated Tap, 5/6 mm
x1



PSANC006
Set Screw Driver
x1



PSANC019
Cortical Awl
x1



/ Revision Screw Driver

PSANC013
Fork Guide
x1



PSANC025
Rod Introducer
x1




PSANC015
Distraction Pliers
x1



PSANC014
Compression Pliers
x1



PSANC016
Rod Caliper
x1



PSANC024
Open Screw Extender
x1



PSANC022
Open Set Screw Guide
x1




PSANC007
Universal Handle
x2



PSANC008
Torque Limiting Handle
x1



PSANC020
Navigation Clamp
x1



Contact Safe Orthopaedics

For more information about our products and company, please visit www.safeortho.com or contact us today.

Western Europe / International

Safe Orthopaedics, SAS
Parc des Bellevues
Allée R. Luxembourg, Le Californie
95610 Eragny sur Oise
FRANCE
Ph: +33 1 34 21 12 43
Fax: +33 1 34 21 12 00

North America

Safe Orthopaedics, LLC
3175 Lenox Park Blvd.
Suite 108
Memphis, TN 38115
USA
Toll Free Ph: (855) 837-3959
Toll Free Fax: (855) 837-3958



SteriSpine™ PS

Pedicle Screw System



Multi-Axial cannulated pedicle screws



**PSANC002
Screw Extender
Ring
x2**



**PSANC005
Screw Driver
x2**



**PSANC017
Retraction
Clamp
x1**



**PSANC001
Screw Extender
Arm
x2**



**PSANC003
Rod Pusher
/ Set Screw Holdex2**



- KITMC530
- KITMC535
- KITMC540
- KITMC545
- KITMC550
- KITMC635
- KITMC640
- KITMC645
- KITMC650
- KITMC655
- KITMC735
- KITMC740
- KITMC745

- 5.5 x 30mm
- 5.5 x 35mm
- 5.5 x 40mm
- 5.5 x 45mm
- 5.5 x 50mm
- 6.5 x 35mm
- 6.5 x 40mm
- 6.5 x 45mm
- 6.5 x 50mm
- 6.5 x 55mm
- 7.5 x 35mm
- 7.5 x 40mm
- 7.5 x 45mm

- Cannulated Screw Kit
- Cannulated Screw Kit
- Cannulated Screw Kit
- Cannulated Screw Kit
- Cannulated Screw Kit
- Cannulated Screw Kit
- Cannulated Screw Kit
- Cannulated Screw Kit
- Cannulated Screw Kit
- Cannulated Screw Kit
- Cannulated Screw Kit
- Cannulated Screw Kit
- Cannulated Screw Kit



Pre-Bent and Straight Rod kits



- KITRMI535
- KITRMI540
- KITRMI550
- KITRMI560
- KITRMI570
- KITRMI580
- KITRMI120
- KITRMI180

- 5.5mm x 35mm
- 5.5mm x 40mm
- 5.5mm x 50mm
- 5.5mm x 60mm
- 5.5mm x 70mm
- 5.5mm x 80mm
- 5.5mm x 120mm
- 5.5mm x 180mm

- Pre-Bent Rod
- Pre-Bent Rod
- Pre-Bent Rod
- Pre-Bent Rod
- Pre-Bent Rod
- Pre-Bent Rod
- Straight Rod
- Straight Rod

Contact Safe Orthopaedics

For more information about our products and company, please visit www.safeortho.com or contact us today.

Western Europe / International

Safe Orthopaedics, SAS
Parc des Bellevues
Allée R. Luxembourg, Le Californie
95610 Eragny sur Oise
FRANCE
Ph: +33 1 34 21 12 43
Fax: +33 1 34 21 12 00

North America

Safe Orthopaedics, LLC
3175 Lenox Park Blvd.
Suite 108
Memphis, TN 38115
USA
Toll Free Ph: (855) 837-3959
Toll Free Fax: (855) 837-3958